



NEW GENERATION ↗

↘ To make the magic happen with your face and body!

To make the magic happen with
your face and body!

Cutting-Edge Technology

+1.3_{Times}
Effectiveness
improved*

+80%
Safety and comfort
improved*

-1.0_{Times}
Treatment time
reduced*


Clinical trials tested on 35 peoples.
*Determined by expert assessment

MICRO-PULSED TECHNOLOGY

Various cartridges up to 12 types

Customized treatment depending on the skin type, treatment area and depth of dermis to subcutaneous fat layer.


FL HP×7types



SD HP×1types




PC HP×4types



Ultra-Fast MP Mode

MPT (Micro-Pulsed Technology) is the latest generation aesthetic solution, combining advanced micro and macro ultrasound energy. Evolved from user feedback, it enables faster treatments with multiple application options.



Micro Pulse Mode




Normal Mode



FLHP

- ◎ 7 types: Face 4 types (1.5mm+2.0mm+3.0mm+4.5mm), Body 3 types (6.0mm+9.0mm+13mm).
- ◎ Each cartridge supports 2 modes (normal & Micro pulsed).
Enable to change mode with a single touch of the screen.



VARIOUS CARE-DEPTH CARTRIDGES

MFU (Microfocused Ultrasound) waves stimulate or ablate tissue through Neocollagenesis process. Similar to a magnifying glass, the focusing mechanism concentrates ultrasound waves at a specific point, resulting in visible outcomes including skin lifting, tightening, and body contouring.

Facial Lifting & Skin Tightening


				
Focal Depth	1.5mm	2.0mm	3.0mm	4.5mm
Mode	MP/Normal	MP/Normal	MP/Normal	MP/Normal

Body Contouring

			
Focal Depth	6.0mm	9.0mm	13.0mm
Mode	MP/Normal	MP/Normal	MP/Normal

SD HP

- ◎ 1 type: Face (1.5mm+2.0mm+3.0mm+4.5mm, can adjust the penetration depth manually).
- ◎ Maximize the clinical results, minimize the treatment pain and swelling. Double lifting effects with SD overlapping treatment after linear and dotting procedure.




“ The physician can adjust the penetration depth manually with a single touch of the screen. ”

PC HP

- ◎ 4 types: Face(1.5mm+2.0mm+3.0mm+4.5mm).
- ◎ Each cartridge supports 2 modes (normal & Micro pulsed). Enable to change mode with a single touch of the screen.

- ▶ Endless operation of the circular
- ▶ Multiple cartridges with different depths of care
- ▶ Multi-mode combined nursing



				
Focal Depth	1.5mm	2.0mm	3.0mm	4.5mm
Mode	MP/Normal	MP/Normal	MP/Normal	MP/Normal

Non-invasive
Facelift, Skin Tightening,
Body Contouring



Wrinkle
Removal



Body
Contouring



Face Lifting
& Tightening



Enhance Jawline
Definition

6 STEPS SKIN TIGHTENING

STEP1

Dermal collagen synthesis decreases progressively over time.

STEP2

Assess the treatment area and design a shot protocol based on individual anatomy and severity of laxity.

STEP3

Apply ultrasound coupling gel to ensure optimal energy transmission, accounting for specific skin conditions.

STEP4

Precisely deliver MFU energy to targeted depths of 1.5mm, 3.0mm, 4.5mm, and 13.0mm for comprehensive tissue coverage.

STEP5

The MFU energy stimulates a controlled wound healing response, initiating phases of coagulation, wound contraction, and subsequent neocollagenesis.


STEP6

This initiates the skin's natural repair and remodeling processes, resulting in improved structural elasticity and firmness.


REAL CASE STUDY RESULTS SPEAKING CLEARLY

The data comes from real and rigorous clinical reports.

Before and After



*1 Treatment MP & Normal Mode



*1 Treatment MP & Normal Mode



*2 Treatment MP & Normal Mode



*2 Treatment MP & Normal Mode